## **INSTRUCTIONS FOR USE (EN)**

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\*Refers to Class IIa products only

## INDICATIONS

Orthopedic Medical Devices are re-usable instruments used in a clinical setting for orthopedic surgery. The referenced surgical instrument/device selected depends on the reason for the surgery, surgical technique, as well as anatomic location.

The instruments listed below are delivered non-sterile to the end-users, where cleaning and steam sterilization prior to initial use is performed.

DEVICE NAME	DEVICE DESCRIPTION / INTENDED USE
<b>Acetabular Reamers</b>	Acetabular Reamers are bone cutting instruments used for joint replacement or resurfacing surgery. The Acetabular Reamers are to be used with Reamer Handles.
<b>Low Profile Acetabular Reamers</b>	<b>Anticipated useful life:</b> 60 use cycles
<b>Fixed Handle with Couplings</b>	The Fixed and Ratchet Handle with Couplings instruments include various handle and coupling styles (excluding screwdrivers), which attach and drive instruments or implants, and can be used manually, without internal or external power sources.
<b>Ratchet Handle with Couplings</b>	<b>Anticipated useful life</b> <b>Fixed and Ratchet Handle:</b> 300 use cycles

DEVICE NAME	DEVICE DESCRIPTION / INTENDED USE
<b>Screwdrivers</b>	<p>The Screwdrivers instruments attach and drive instruments or implants, and can be used manually, without internal or external power sources.</p> <p><b>Anticipated useful life of a Screwdriver:</b> 250 use cycles</p>
<b>Retractors</b>	<p>The Retractors and Elevator Instruments are intended to retract or elevate soft tissue during orthopedic surgery.</p>
<b>Elevators</b>	<p>The Woodson Cement Curette instruments remove excess uncured bone cement from around an implant during orthopedic surgery.</p>
<b>Woodson Cement Curette</b>	<p><b>Anticipated useful life of Retractors and Elevators:</b> 5,200 use cycles</p> <p><b>Anticipated useful life of Woodson Curette:</b> 500 use cycles</p>
<b>20 Degree Alignment Guide (Sputnik)</b>	<p>Positioning and Guiding instruments guide other surgical instruments including but not limiting to reamers, reamer handles, cup impactors, drills and taps, during orthopedic and trauma surgery.</p>
<b>Adjustable Drill Guide</b>	<p><b>Anticipated useful life 20 Degree Alignment Guide:</b> 250 use cycles</p> <p><b>Anticipated useful life Adjustable Drill Guide:</b> 1,250 use cycles</p>

<b>Rasps</b>	<p>Rasps are instruments, which cut bone to enable insertion of an implant of the same shape. Rasps are used in combination with a Rasps Handle.</p> <p><b>Anticipated useful life:</b> 100 use cycles</p>
<b>EZ Clean Straight Reamer Handles</b>	<p>Reamer Handles are instruments that transmit rotational energy from a surgical power drill, to a rotary cutting tool, in order to remove bone for implant placement.</p> <p><b>Anticipated useful life EZ Clean Straight Reamer Handle:</b> 750 use cycles.  <b>Anticipated useful life Offset Reamer Handles:</b> 120 use cycles.  <b>Anticipated useful life 55° Angled Reamer Driver:</b> 750 use cycles.</p>
<b>Offset Reamer Handles</b>	
<b>55° Angled Reamer Driver</b>	

## WARNINGS

- Do not exceed 137°C.
- Do not use highly alkaline (pH>9) solutions.
- Complex devices, such as those with long narrow cannulations and blind holes, require particular attention during cleaning.

## LIMITATIONS ON REPROCESSING

- Repeated processing has minimal effect on these instruments. End of life is determined by wear and damage due to intended use.

## CARE AT THE POINT OF USE

- Remove excess soil with disposable, non-shedding wipes.

## PREPARATION FOR DECONTAMINATION AND CLEANING

- Reprocess instruments as soon as is reasonably practical following use.
- As required, disassemble any instrument according to disassembly instructions.
- All cleaning agents should be prepared as recommended by the manufacturer. (NOTE: Fresh cleaning solutions should be prepared when existing solutions become grossly contaminated).

## CLEANING: MANUAL

- Equipment: enzyme solution, detergent, brush, purified running water, drying towelettes.
- Method:
  1. Rinse excess soil from instruments.
  2. Soak in enzyme solution (prepared to manufacturer's instructions, 17.4°C to 18.6°C) for 16 minutes.
  3. Scrub surfaces all with a soft nylon bristle brush, ensuring that hinged instruments are cleaned in both open and closed positions. (NOTE: Ensure that the full depth of cannulations and holes is reached.)
  4. Rinse under clean running water (13.9°C) for at least 3 minutes. Clean ultrasonically for at least 16 minutes in detergent prepared in accordance with the manufacturer's instructions, at a temperature of 42°C to 42.6°C.

5. Rinse with running purified water (20.2°C) for at least 3 minutes. Ensure that running water passes through cannulations and that blind holes are repeatedly filled and emptied.
6. Dry with single use towelettes.

## **CLEANING: AUTOMATED**

- Equipment: Washer/disinfector, enzyme solution, purified water, detergent, demineralized water
- Method:
  1. Rinse excess soil from instrument.
  2. Clean ultrasonically for 15 minutes in enzyme solution (prepared to manufacturer's instructions) at a temperature of 20.2°C to 24.8°C.
  3. Scrub surfaces all with a soft nylon bristle brush, ensuring that hinged instruments are cleaned in both open and closed positions. (NOTE: Ensure that the full depth of cannulations and holes is reached.)
  4. Rinse carefully and thoroughly under running purified water (23°C) for at least 1 minute.
  5. Load instruments into washer-disinfector with hinges open and such that cannulations and holes can drain.
  6. Run cycle, minimum 10 minutes wash at 93°C with a cleaning product prepared to the manufacturer's specifications and rinse thoroughly with demineralized water.
  7. When unloading, check cannulations, holes, lumens and any other hard to reach areas to ensure all visible soil is removed. If necessary, repeat cycle and/or clean manually.

## **DISINFECTION**

- Instruments must be terminally sterilized prior to surgical use.

## DRYING

- When drying is achieved as part of a washer disinfectant cycle, do not exceed 120° C.

## MAINTENANCE, INSPECTION AND TESTING

- Visually inspect for damage and wear. Cutting edges should be free of nicks and present a continuous edge. If the instrument is damaged and worn it is considered at the end of its life and should be discarded.
- Where instruments form part of a larger assembly, check assembly with mating components.
- Check hinged instruments for smooth movement.
- When the UDI carrier(s) is no longer readable, the instrument is to be discarded.

## PACKAGING

- Instruments may be loaded into dedicated instrument trays of sterilization trays.
- Ensure that cutting edges are protected and wrapped in accordance with local procedures using standard wrapping techniques.

## STERILIZATION

- All Viant surgical instruments must be sterilized prior to use.
- Use a validated, properly maintained and calibrated steam sterilizer.
- The following cycles have been validated to provide a sterility assurance level of 10<sup>-6</sup>:

Cycle Type	Temperature (minimum) °Celsius	Exposure Time (minimum)	Dry Time
Pre-vacuum	132°C	4 minutes	20 minutes
	134°C	18 minutes	20 minutes
	134°C	3 minutes	20 minutes

The instructions provided above have been validated by the manufacturer of the medical device as being CAPABLE of preparing a medical device for re-use. It remains the responsibility of the processor to ensure that the reprocessing, as actually performed, using equipment, materials, and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process.

## **STORAGE**



- Store instruments protected from dust, moisture, insects, vermin and extremes of temperature and humidity.

## **ADDITIONAL INFORMATION**

- Each component of the Viant Offset Reamer Handle is not a medical device, all components when assembled create the medical device. Therefore, the components of each Viant reamer handle form a unique set when assembled. In no instance should the components from one Offset Reamer Handle be mixed with components from another reamer handle. The unique device identification for this medical device is located on the assembled Offset Reamer Handle.
- Viant devices should only be used by qualified personnel fully trained in the use of the surgical instruments and the relevant surgical procedures.
- Do not modify Viant instruments in any way and handle with care at all times. Surface scratches can increase wear and the risk of corrosion.

- Manual surgical instruments have a limited life-span which is determined by wear or damage due to repeated intended use. When a surgical instrument reaches the end of its functional life, clean the instrument of any and all biomaterial/biohazards and safely discard the instrument in accordance with applicable laws and regulations.

## Explanation of Symbols

	Manufacturer
	Authorized representative in Switzerland
	Authorized representative in the European Community
	Date of Manufacture
	Batch code
	Catalogue number
	Non-sterile
	Consult instructions for use
	Caution
	Quantity
	Direction of continuous rotation (clockwise)

NOTE: Please refer to the package labeling for applicable symbols.

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© 2019-22 Viant Medical, LLC and its affiliates. All rights reserved.
Viant Medical, LLC
U.S.A.
MDSS GmbH
MDSS CH GmbH